



April 15, 2021

Terahertz Nanoscanning Spectrometer & 3D Imager (TNS3DI) **A new tool for semiconductor metrology**

The Applied Research & Photonics model TNS3DI – Terahertz Nanoscanning Spectrometer and 3-D Imaging system, utilizing its proprietary Dendrimer Dipole Excitation technology to deliver what we have demonstrated is unsurpassed analytical, and imaging capability to enable the revolution in advanced materials and structures. The system has demonstrated the ability to pinpoint and catalyze the deep investigation and understanding scientists and engineers all over the world have been searching for, to solve some of the most challenging problems in advanced materials, semiconductors, and related devices, photonics and optoelectronics including displays, medical, biological and pharmaceutical

applications, and atomic/nuclear level investigations. Two models are available: (1) TNS3DI-FC, fiber-coupled beam delivery for the most versatility, X, Y = 200 mm standard, up to 450 mm option. (2) TNS3DI-Fs, free-space beam delivery for small samples, X, Y ≤ 100 mm. The TNS3DI delivers combined capabilities of AFM, SEM, TEM, and Time-domain Spectrometer. It is also a diagnostic tool to be used in for the US Military and government laboratories who rely heavily on semiconductors (microchips) in day-to-day operations.





Semiconductor and nanomaterials

Terahertz radiation (T-ray) can penetrate semiconductors and most non-metals, so, one can probe both the surface and sub-surface regions in a non-destructive and non-contact route. Layer-by-layer imaging allows inspection of layers on a wafer one at a time. A few examples cited below showing the strengths. Other test protocols may be developed.




Fundamental breakthrough

ARP’s Terahertz Scanning Spectrometer (TeraSpectra) is a Terahertz Nano-Scanning Spectrometer/3D Imaging system that has a few key technology innovations:


- 1) A new terahertz generation mechanism, dendrimer dipole excitation (DDE), for high power, up to 30 THz, continuous wave, stable terahertz generation.
- 2) Overcoming the Abbe diffraction limit for lattice resolution image generation with bigger (terahertz) wavelength.
- 3) Replacing many functionalities of AFM/SEM/TEM by T-ray technique, it uniquely identifies location, depth, and type of defects, where it exists.
- 4) Only technology available to see interior (sub-surfaces) in a non-destructive route with layer-by-layer imaging and analysis

Currently, there is no measurement technology that has the capability to provide an equivalent richness of information that ARP’s *TNS3DI* system can deliver *without* damage or destruction of the test sample. A few examples of specific applications are given below. Many more will be invented by the end users and experts in their respective areas.

| Applications for the semiconductor and nanotechnology industry | |
|--|--|
| Test specification | Exemplary Results |
| <p>Whole wafer inspection Whole wafer may be imaged at high speed and each die may be inspected separately Terahertz image of unpatterned metalized area shows higher reflection (bluish).</p> | <p>The figure shows two 3D surface plots of a wafer. The left plot shows a regular grid of peaks, with a green arrow pointing to a peak. The right plot shows a similar grid but with a significant dip in the center, also with a green arrow pointing to it. Below each plot is a line graph showing the 'Grey Value' vs 'Distance (µm)'. The right graph is labeled 'Good' and 'bad'.</p> |


| | | |
|--|--|--|
| <p>Lattice resolution imaging Lattice image of metallic nickel shown. Lattice parameter was measured: 0.353 nm.</p> |  | |
| <p>3D (volumetric) imaging Here a 3D image of alumina shown over (50 nm)³ volume</p> |  | |
| <p>Phase contrast imaging via short-time Fourier transform frequency spectrum</p> | <p style="text-align: center;">H2@C60 Short-Time Fourier Transform Frequency Spectrum</p>  | |

SiC wafer's surface inspection. T-ray sees more than anything else.



Top left: Optical image of a SiC epitaxial 100 mm whole wafer. Top right: T-ray image of the same wafer
 Bottom: two different views of a zoomed in segment. Many features are seen from T-ray image that cannot be seen otherwise.


Layer by layer analysis via 3D imaging.
 A volume image may be sliced in to any number of layers on any of the orthogonal planes.
 XY-plane and XZ-planes are shown



Semiconductor interfaces. Deposition involves transforming a solid material in to vapor phase and then back in to solid by the deposition process. The deposited materials' lattice will/may suffer from defects such as stacking fault and dislocations. Volume and 2D image across the interface allow critical metrology of parameters such as layer thickness and lattice defects.


Different kind of interfaces

- Misfit-lattice interface (type-1)
- Rough surface interface (type-3) → Layer-1/Layer-2
- Diffusion type interface (type-4) → Layer-2/Layer-3




A single face 2D image across the depth, 200 nm × 200 nm. Sharp interface line runs across between the layers. The layer thicknesses may be measured.

- Identical lattice interface (type-2)




Close-up 3D image of a gold nanoparticle



Left: Close-up of a gold nanoparticle lattice image (10 nm)².
 Right: Further close-up of the same over (1 nm)².

© Anis Rahman, PhD, Applied Research & Photonics for T-ray Image a Day


Very high resolution image of a gold nanoparticle



Left: Very high resolution image of a gold nanoparticle
 Right: A different view of the same.

© Anis Rahman, PhD, Applied Research & Photonics for T-ray Image a Day


Size and distribution of any nanoparticle in dry or wet form.
 Measured quantum dot size ~7.7 nm





| | |
|---|--|
| <p>Non-destructive, nano-scale testing and metrology down to <1 nm feature on and under the surface of a wafer.</p> | |
| <p>Size and distribution of any nanoclusters in dry or wet form. Silica nanoparticles in isopropyl alcohol (nanoslurry) showing the distribution of nanoclusters in suspension. Nanoparticle size is ~10 nm Nanocluster size is ~17.5 nm</p> | |
| <p>Measurement of the nanovoids Measured nanovoid size ~2.8 nm</p> | |

| | |
|--|--|
| <p>Epitaxial semiconductor layer thickness</p> <p>Measured SiGe layer thickness ~18 nm</p> | <p>The figure shows a cross-sectional HRTEM image of a SiGe layer. The image displays wavy lattice fringes. A line scan graph to the right plots Gray Value (0 to 150) against Distance (nm) (0 to 30). A sharp dip in the gray value is observed at approximately 18 nm, corresponding to the SiGe layer thickness.</p> |
| <p>Dislocation in epitaxial semiconductor</p> <p>Dislocations are a type of defect in crystalline structures where the atoms are out of position in the lattice.</p> <p>Dislocations are generated and moved as a result of an applied stress.</p> | <p>The figure shows an HRTEM image of a dislocation in an epitaxial semiconductor. A dislocation core is visible as a localized distortion in the lattice fringes. A line scan graph to the right shows the corresponding intensity profile, with a sharp peak at the dislocation core.</p> |
| <p>Lattice damage due to electric stress.</p> <p>Here a GaN grown on silicon damaged by high electric field stress</p> | <p>The figure shows an HRTEM image of GaN grown on silicon. The image displays a clear lattice structure. A region of the lattice is labeled "Damaged lattice area" and shows distorted fringes, while the rest of the image is labeled "Good lattice area" and shows regular fringes.</p> |

Graphene exfoliates' 3D imaging
 Number of layers in an exfoliate is measured
 Thickness of each layer is quantified




Roughness of metalized wafer.
 Optical image does not show the roughness




T-ray image of the same wafer on the left. Defects are visible. Surface roughness may be quantified.

Surface roughness quantification




Measurement of surface roughness from T-ray image, $\pm 0.25 \mu\text{m}$

Deep-level terahertz spectroscopy for detection of minute changes and/or molecular identification.
 $t = 0 \text{ nm} \rightarrow$ surface, Ge
 $t = 600 \text{ nm} \rightarrow$ Si <100>
 Shifts in peaks for Ge wrt Si
 Peaks at 0.6 THz and 1.3 THz are due to ambient moisture.




Resistivity map of the whole wafer – may be extended up to 450 mm



Resistivity map (contour plot) of 100 mm wafer

Force scan (friction measurement). Demonstrated with a patterned semiconductor wafer and human stratum corneum.



Top left: Image of a segment of a 200 mm patterned wafer. Bottom left: the frictional force variations as a function of scan distance.

Top right: Human stratum corneum treated with FXII protein
 Bottom right: A lateral force scan (friction measurements).

Contact info

Anis Rahman, PhD, Phone: +1 (717) 623-8201. Email: a.rahman@arphotronics.net
 Applied Research & Photonics, Inc., 470 Friendship Road, Suite 10, Harrisburg, PA 17111,
<http://arphotronics.net>